

BioSpherix

XVIVO SYSTEM

Closed Cell Incubation and Processing

FIRST AND ONLY BARRIER ISOLATOR OPTIMIZED FOR CELLS



Cells are totally isolated in closed system. No bunny suits required.



Incubators and other equipment open only into optimized and aseptic workspace.

XVIVO SYSTEM - a practical and economical alternative to cleanrooms for cGMP compliant production of human cells for clinical use.

Better Process Control

- Full time optimization of Temperature, RH, CO₂/pH, O₂
- More consistent phenotype

Better Contamination Control

- Closed protects cells from people and room
- Closed protects people from cells
- Isolation protects Sample A from Sample B

Better Quality Control

- PAT/QbD parametric release
- PAT/QbD continuous process improvement

Better Pocketbook Control

- 10% of total cost of ownership
- No brick-&-mortar risks (flexible & movable)